





User Guide: Getting Started



Please select **“Sign up now”** to create a DREX account

(3)

It is important to choose a username that is unrelated to yourself, for anonymity purposes



Enter your details.
Submit by clicking **“Register”**



(4)

The screenshot shows a 'Terms of Use' form. At the top, it reads: 'THIS AGREEMENT IS MADE BETWEEN THE END USER OF THIS SOFTWARE AND THE UNIVERSITY OF DURHAM, WHOSE REGISTERED ADDRESS IS THE PALATINE CENTRE, STOCKTON ROAD, DURHAM, DH1 3LE.' Below this, it says 'WHEREAS: A. Durham Reading and Exploration training (DREX) is self-adjusting, software designed by researchers at the University of Durham (hereinafter "the University") with the aim of attempting to assist in the rehabilitation of individuals who suffer from homonymous visual field'. At the bottom, there is a paragraph: 'Your name is used for registration purposes only. All personal data and performance data will be stored anonymously / not identifiable. You may optionally choose to share your personal and performance data with a doctor of your choice, and revoke such access at any time.' At the very bottom are two buttons: 'Accept' and 'Decline'. A speaker icon is circled in red, with a dotted line pointing to it from the right.

For assistance the text may be read aloud by clicking the speaker icon

Read and agree to terms and conditions of use

(5)

The screenshot shows a 'Details of your condition' form. It starts with the instruction: 'Please provide some details of your condition.' The first question is 'Which side of your visual field is affected?' with radio button options for 'Left', 'Right', and 'Both'. The 'Right' option is selected. Below this is a section for 'Your gender' with radio button options for 'Male' and 'Female'. The word 'optional' is written next to the gender section. At the bottom, there is a text input field with the label 'Describe the cause behind your current condition' and the word 'optional' next to it. The question text and the 'optional' label for the gender section are circled in red.


This field is required, to enable DREX to benefit your rehabilitation

Other details are optional; you do not need to provide further information to use DREX



User Guide: Initial Assessment

1a) Visual Search Tasks




Identify and click on the pen in a scene of every day items

1b) Visual Search Tasks

Identify and click the numbers 1-20 in sequence

12 18 7 4
10 6 2 20
17 3 19 5
8 1 13 15
14 9 11 16


2) Short Term Memory Task



You will presented with a set of numbers for 10 seconds.
Then a second set of numbers. Click the number that appeared in the first set.



3) Quality of Life Task



You will be asked to rate how difficult you find specific tasks. **Slide the green bar to respond**

When satisfied with response click **“Next Question”**

You will be asked to confirm your response is correct. Click **“Yes”** if you are happy to proceed.


4) Reading Task

The man stood looking on for a while. He began to be tired, and hungry too, for he had given away his last penny in his joy at getting the cow. At last he could go no further, for the stone tired him sadly: and he dragged himself to the side of a river, so he could take a drink of water and rest a while. So he laid the stone carefully by his side on the bank: but, as he stooped down to drink, he forgot it, pushed it a little, and down it rolled, plump into the stream.

Read the passage presented. Click the screen to continue.

You will be asked 3 questions about the passage you have just read.

5) Perimetry Task




Looking at the red dots in the centre of the screen, click the location where a different colour dot appears. Tap the space bar for the next screen if no dot appears.



User Guide: Training

Exploration



These exercises require you to identify the odd one out on the page. There are 3 types of task; to identify and CLICK on the different colour, size or shape.

Reading

coffee <i>Word</i> <i>Tap the space</i>	fefcoe <i>Non-Word</i> <i>Click the screen</i>	doffee <i>Non-Word</i> <i>Click the screen</i>
--	---	---

This exercise requires you to identify if the image shows a word or a non-word. For words tap the SPACE bar, for non-words click the screen.